



## دوربین تحقیقاتی میکروسکوپ

### Microbin 2



Model	Microbin 2
Sensor Type	Sony CMOS
Sensor Size	1/2.9"
Resolution	2MP, 1920x1080
Pixel Size	2.8μm x 2.8μm
Frame	30fps
Color/Mono	Color
Shutter Mode	Rolling
Exposure Mode	Auto/Manual
Picture Format	TIFF/JPG/PNG/DICOM
Operating System	Windows 7/8/10(32/64bit),Mac
PC Software	Mosaic V2
SDK	C/C++, C#, Directshow
Optical Interface	C Mount
Data Interface	USB2.0
Operating Environment	Temperature: 0-40°C Humidity: 10%-85%
Camera Size	68X68X42.5(mm)
Camera Weight	236g