


clear & amber
modular burettes
 class "AS"

- All components of the modular burette are replaceable.
- Multiple burette tubings can be used with one single stopcock.
- Modular construction with detachable burette tube and stopcock enables fast dismantle and easy cleaning.
- Inscriptions and graduations are printed with special **blue enamel**.
- Clear burettes are manufactured with schellbach stripe
- Other technical specifications are exactly the same as other types of burettes.

Modular construction with detachable burette tube and stopcock enables fast dismantle and easy cleaning.



clear - schellbach blue graduations	amber - w/o schellbach white graduations	volume	tolerance	subdivision	quantity in box
017.17.010	017.24.010	10 ml	± 0.020 ml	0.02 ml	1 piece
017.17.025	017.24.025	25 ml	± 0.030 ml	0.05 ml	1 piece
017.17.050	017.24.050	50 ml	± 0.050 ml	0.10 ml	1 piece

burette tubes for modular burettes

clear - schellbach blue graduations	amber - w/o schellbach white graduations	volume	tolerance	subdivision	quantity in box
017.25.010	017.31.010	10 ml	± 0.020 ml	0.02 ml	1 piece
017.25.025	017.31.025	25 ml	± 0.030 ml	0.05 ml	1 piece
017.25.050	017.31.050	50 ml	± 0.050 ml	0.10 ml	1 piece